

### Proteins

5 The present invention relates to proteins derived from *Streptococcus agalactiae*, nucleic acid molecules encoding such proteins, and the use of the proteins as antigens and/or immunogens and in detection/diagnosis. It also relates to a method for the rapid screening of bacterial genomes to isolate and characterise bacterial cell envelope associated or secreted proteins.

10 The Group B Streptococcus (GBS) (*Streptococcus agalactiae*) is an encapsulated bacterium which emerged in the 1970s as a major pathogen of humans causing sepsis and meningitis in neonates as well as adults. The incidence of early onset neonatal infection during the first 5 days of life varies from 0.7 to 3.7 per 1000 live births and causes mortality in about 20% of cases. Between 25-50% of neonates surviving early onset infections frequently suffer neurological sequelae. Late onset neonatal  
15 infections occur from 6 days to three months of age at a rate of about 0.5 - 1.0 per 1000 live births.

20 There is an established association between the colonisation of the maternal genital tract by GBS at the time of birth and the risk of neonatal sepsis. In humans it has been established that the rectum may act as a reservoir for GBS. Susceptibility in the neonate is correlated with the a low concentration or absence of IgG antibodies to the capsular polysaccharides found on GBS causing human disease. In the USA strains isolated from clinical cases usually belong to capsular serotypes Ia, Ib, II, III although serotype V may be of increasing significance. Type VIII GBS is the major  
25 cause of neonatal sepsis in Japan.

A possible means of prevention involves intra or postpartum administration of antibiotics to the mother but there are concerns that this might lead to the emergence

of resistant organisms and in some cases allergic reactions. Vaccination of the adolescent females to induce long lasting maternally derived immunity is one of the most promising approaches to prevent GBS infections in neonates. The capsular polysaccharide antigens of these organisms have attracted most attention as with  
5 regard to vaccine development. Studies in healthy adult volunteers have shown that serotype Ia, II and III polysaccharides are non-toxic and immunogenic in approximately 65%, 95% and 70% of non-immune adults respectively. One of the problems with using capsule antigens as vaccines is that the response rates vary according to pre-immunisation status and the polysaccharide antigen and not all  
10 vaccinees produce adequate levels of IgG antibody as indicated in vaccination studies with GBS polysaccharides in human volunteers.

Some people do not respond despite repeated stimuli. These properties are due to the T-independent nature of polysaccharide antigens. One strategy to enhance the  
15 immunogenicity of these vaccines is to enhance the T cell dependent properties of polysaccharides by conjugating them to a protein. The use of polysaccharide conjugates looks promising but there are still unresolved questions concerning the nature of the carrier protein. A conjugate vaccine against GBS would require at least 4 different conjugates to be prepared adding to the cost of a vaccine.

20 Approaches to vaccination against GBS infections which rely on the use of capsular polysaccharides have the disadvantage that response rates are likely to vary considerably according to pre-immunisation status and the particular type of polysaccharide antigen used. Results of trials with conjugate vaccines in human  
25 volunteers have indicated that response rates may only be around 65% for some of the key capsule antigens (Larsson *et al.*, *Infection and Immunity* 64:3518-3523 (1996)). It is also not clear whether all individuals responding to the vaccine would have adequate levels of polysaccharide specific IgG which can cross the placenta and

afford immunity to neonates. By conjugating a protein carrier to the polysaccharide antigen it may be possible to convert them to T-cell dependent antigens and enhance their immunogenicity.

- 5 Preliminary studies with GBS type III polysaccharide-tetanus toxoid conjugate have been encouraging (Baker *et al.*, *Reviews of Infectious Diseases* 7:458-467 (1985), Baker *et al.*, *The New England Journal of Medicine* 319:1180-1185 (1988), Paoletti *et al.*, *Infection and Immunity* 64:677-679 (1996), Paoletti *et al.*, *Infection and Immunity* 62:3236-3243 (1994)) but in developed countries the use of tetanus may be  
10 disadvantageous since most adults will have been immunised against tetanus within the past five years. Additional boosters with tetanus toxoid may cause adverse reactions (Boyer., *Current Opinions in Pediatrics* 7:13-18 (1995)). The polysaccharide conjugate vaccines have the disadvantage of being costly to produce and manufacture in comparison with many other kinds of vaccines. There is also the  
15 possible risk of problems caused by the cross reactivity between GBS polysaccharides and sialic acid-containing human glycoproteins.

- Recent evidence suggests that bacterial surface proteins also may be useful to confer immunity. A protein called Rib which is found on most serotype III strains but rarely  
20 on serotypes Ia, Ib or II confers immunity to challenge with Rib expressing GBS in animal models (Stalhammar-Carlemalm *et al.*, *Journal of Experimental Medicine* 177:1593-1603 (1993)). Another surface protein of interest as a component of a vaccine is the alpha antigen of the C proteins which protected vaccinated mice against lethal infection with strains expressing alpha protein. The amount of this  
25 antigen expressed by GBS strains varies markedly, however an alternative to polysaccharides as antigens is the use of protein antigens derived from GBS. Recent evidence suggest that the GBS surface associated proteins Rib and alpha C protein may be used to confer immunity to GBS infections in experimental model systems

(Stalhammar-Carlemalm *et al.*, (1993) [*supra*], Larsson *et al.*, (1996) [*supra*]). However these two proteins are not conserved in all serotypes of GBS which cause disease in humans. Assuming that these antigens would be immunogenic and elicit protective level responses in humans they would not confer protection against all  
5 infections caused by GBS as 10% of infectious Group B streptococci do not express Rib or C protein alpha.

This invention seeks to overcome the problem of vaccination against GBS by using a novel screening method specifically designed to identify those Group B  
10 Streptococcus genes encoding bacterial cell surface associated or secreted proteins. The proteins expressed by these genes may be immunogenic, and therefore may be useful in the prevention and treatment of Group B Streptococcus infection. For the purposes of this application, the term immunogenic means that these proteins will elicit a protective immune response within a subject. Using this novel screening  
15 method a number of genes encoding novel Group B Streptococcus proteins have been identified.

Thus in a first aspect, the present invention provides a Group B Streptococcus protein, polypeptide or peptide having a sequence selected from those shown in  
20 figure 1, or fragments or derivatives thereof.

It will be apparent to the skilled person that proteins and polypeptides included within this group may be cell surface receptors, adhesion molecules, transport proteins, membrane structural proteins, and/or signalling molecules.

25 Alterations in the amino acid sequence of a protein can occur which do not affect the function of a protein. These include amino acid deletions, insertions and substitutions and can result from alternative splicing and/or the presence of multiple translation

5 Thus, the present invention includes derivatives or variants of the proteins, polypeptides, and peptides of the present invention which show at least 50% identity to the proteins, polypeptides and peptides described herein. Preferably the degree of sequence identity is at least 60% and preferably it is above 75%. More preferably still it is above 80%, 90% or even 95%.

The term identity can be used to describe the similarity between two polypeptide sequences. A software package well known in the art for carrying out this procedure is the CLUSTAL program. It compares the amino acid sequences of two polypeptides and finds the optimal alignment by inserting spaces in either sequence as appropriate. The amino acid identity or similarity (identity plus conservation of amino acid type) for an optimal alignment can also be calculated using a software package such as BLASTx. This program aligns the largest stretch of similar sequence and assigns a value to the fit. For any one pattern comparison several regions of similarity may be found, each having a different score. One skilled in the art will appreciate that two polypeptides of different lengths may be compared over the entire length of the longer fragment. Alternatively small regions may be compared. Normally sequences of the same length are compared for a useful comparison to be made.

25 Manipulation of the DNA encoding the protein is a particularly powerful technique for both modifying proteins and for generating large quantities of protein for purification purposes. This may involve the use of PCR techniques to amplify a desired nucleic acid sequence. Thus the sequence data provided herein can be used to

design primers for use in PCR so that a desired sequence can be targeted and then amplified to a high degree.

Typically primers will be at least five nucleotides long and will generally be at least ten  
 5 nucleotides long (e.g. fifteen to twenty-five nucleotides long). In some cases primers  
 of at least thirty or at least thirty-five nucleotides in length may be used.

As a further alternative chemical synthesis may be used. This may be automated.  
 Relatively short sequences may be chemically synthesised and ligated together to  
 10 provide a longer sequence.

Thus in a further aspect, the present invention provides, a nucleic acid molecule  
 comprising or consisting of a sequence which is:

- 15 (i) any of the DNA sequences set out in figure 1 herein or their RNA  
 equivalents;
- (ii) a sequence which is complementary to any of the sequences of (i);
- (iii) a sequence which codes for the same protein or polypeptide, as those  
 sequences of (i) or (ii);
- (iv) a sequence which shows substantial identity with any of those of (i),  
 20 (ii) and (iii); or
- (v) a sequence which codes for a derivative or fragment of a nucleic acid  
 molecule shown in figure 1.

The term identity can also be used to describe the similarity between two individual  
 DNA sequences. The 'bestfit' program (Smith and Waterman, *Advances in applied*  
 25 *Mathematics*, 482-489 (1981)) is one example of a type of computer software used to  
 find the best segment of similarity between two nucleic acid sequences, whilst the  
 GAP program enables sequences to be aligned along their whole length and finds the  
 optimal alignment by inserting spaces in either sequence as appropriate.

The present invention includes nucleic acid sequences which show at least 50% identity to the nucleic acid sequences described herein. Preferably the degree of sequence identity is at least 60% and preferably it is above 75%. More preferably still it is above 80%, 90% or even 95%.

5

The term 'RNA equivalent' when used above indicates that a given RNA molecule has a sequence which is complementary to that of a given DNA molecule, allowing for the fact that in RNA 'U' replaces 'T' in the genetic code. The nucleic acid molecule may be in isolated, recombinant or chemically synthetic form.

10

DNA constructs can readily be generated using methods well known in the art. These techniques are disclosed, for example in J. Sambrook *et al*, *Molecular Cloning 2<sup>nd</sup> Edition*, Cold Spring Harbour Laboratory Press (1989). Modifications of DNA constructs and the proteins expressed such as the addition of promoters, enhancers, signal sequences, leader sequences, translation start and stop signals and DNA stability controlling regions, or the addition of fusion partners may then be facilitated.

15

Normally the DNA construct will be inserted into a vector which may be any suitable vector, including plasmid, virus, bacteriophage, transposon, minichromosome, liposome or mechanical carrier. The expression vectors of the invention are DNA constructs suitable for expressing DNA which encodes the desired protein product which may include: (a) a regulatory element (e.g. a promoter, operator, activator, repressor and/or enhancer), (b) a structural or coding sequence which is transcribed into mRNA and (c) appropriate transcription, translation, initiation and termination sequences. The vector may further comprise a selectable marker, for example antibiotic resistance, which facilitates the selection and/or identification of cells containing the vector.

20

25

Expression of the protein is achieved by the transformation or transfection of the vector into a host cell which may be of eukaryotic or prokaryotic origin. For the production of recombinant protein, expression may be inducible expression or expression only in certain types of cells or both inducible and cell-specific. Particularly preferred among inducible vectors are vectors that can be induced for expression by environmental factors that are easy to manipulate, such as temperature and nutrient additives. A variety of suitable vectors, including constitutive and inducible expression vectors for use in prokaryotic and eukaryotic hosts, are well known and employed routinely by those skilled in the art.

A great variety of expression vectors can be used to express the Group B Streptococcus protein(s) of the invention. Such vectors include, among others, chromosomal, episomal and virus-derived vectors, for example, vectors derived from bacterial plasmids, from bacteriophage, from transposons, from yeast elements, from viruses such as baculoviruses, papova viruses, such as SV40, vaccinia viruses, adenoviruses and retroviruses, and vectors derived from combinations thereof, such as those derived from plasmid and bacteriophage genetic elements, such as cosmids and phagemids, all may be used in accordance with the invention. Generally, any vector suitable to maintain, propagate or express nucleic acid to express a polypeptide in a host may be used for expression in this regard. Such vectors thus form yet a further aspect of the invention.

The appropriate DNA sequence may be inserted into the vector by any of a variety of well-known and routine techniques.

The nucleic acid sequence in the expression vector is operatively linked to appropriate expression control sequence(s) including, for instance, a promoter to



direct mRNA transcription. Representatives of such promoters include, but are not limited to, the phage lambda PL promoter, the T3 and T7 promoters, the *E.coli* lac, trp, tac, and  $\lambda$ P<sub>L</sub> promoters, the microbial eukaryote GAL, glucoamylase and cellobiohydrolase promoters and the mammalian metallothionein (mouse) and heat-shock (human) promoters.

In general, expression vectors will contain sites for transcription initiation and termination, and, in the transcribed region, a ribosome binding site for translation. The coding portion of mature transcripts expressed by the constructs will generally include a translation initiating AUG at the beginning and a termination codon appropriately positioned at the end of the polypeptide to be translated.

Representative examples of appropriate hosts for recombinant expression of the Group B Streptococcus protein(s) of the invention include bacterial cells, such as *streptococci*, *staphylococci*, *E.coli*, *streptomyces* and *Bacillus subtilis* cells; fungal cells, such as yeast cells and *Aspergillus* cells; insect cells such as *Drosophila* S2 and *Spodoptera* Sf9 cells; animal cells such as CHO, COS, HeLa and Bowes melanoma cells; and plant cells. Such host cells form yet a further aspect of the present invention.

Microbial cells employed in the expression of proteins can be disrupted by any convenient method, including freeze-thaw cycling, sonication, mechanical disruption, or use of cell lysing agent, such methods which are known to those skilled in the art.

The polypeptide can be recovered and purified from recombinant cell cultures by well-known methods including ammonium sulphate or ethanol precipitation, acid extraction, anion or cation exchange chromatography, phosphocellulose,

chromatography, hydrophobic interaction chromatography, affinity chromatography, hydroxylapatite chromatography and lectin chromatography. Well known techniques for refolding protein may be employed to regenerate active conformation when the polypeptide is denatured during isolation and or purification.

5

The Group B Streptococcus proteins described herein can additionally be used as target antigens to raise antibodies, or to generate affibodies. These can be used to detect Group B Streptococcus.

10 Thus in a further aspect the present invention provides, an antibody, affibody, or a derivative thereof which binds to any one or more of the proteins, polypeptides, peptides, fragments or derivatives thereof, as described herein.

Antibodies within the scope of the present invention may be monoclonal or polyclonal.

15 Polyclonal antibodies can be raised by stimulating their production in a suitable animal host (e.g. a mouse, rat, guinea pig, rabbit, sheep, goat or monkey) when a protein as described herein, or a homologue, derivative or fragment thereof, is injected into the animal. If desired, an adjuvant may be administered together with the protein. Well-known adjuvants include Freund's adjuvant (complete and incomplete) and aluminium  
20 hydroxide. The antibodies can then be purified by virtue of their binding to a protein as described herein and by many other means well-known to those skilled in the art.

Monoclonal antibodies can be produced from hybridomas. These can be formed by fusing myeloma cells and spleen cells which produce the desired antibody in order to  
25 form an immortal cell line. Thus the well-known Kohler & Milstein technique (*Nature* 256 (1975)) or subsequent variations upon this technique can be used.

Techniques for producing monoclonal and polyclonal antibodies that bind to a particular polypeptide/protein are now well developed in the art. They are discussed in standard immunology textbooks, for example in Roitt *et al*, *Immunology* second edition (1989), Churchill Livingstone, London.

5

In addition to whole antibodies, the present invention includes derivatives thereof which are capable of binding to proteins etc as described herein. Thus the present invention includes antibody fragments and synthetic constructs. Examples of antibody fragments and synthetic constructs are given by Dougall *et al* ., *Tibtech* 12 372-379 (September 10 1994).

15

Antibody fragments include, for example, Fab, F(ab')<sub>2</sub> and Fv fragments. Fv fragments can be modified to produce a synthetic construct known as a single chain Fv (scFv) molecule. This includes a peptide linker covalently joining V<sub>h</sub> and V<sub>l</sub> regions, which contributes to the stability of the molecule. Other synthetic constructs that can be used include CDR peptides. These are synthetic peptides comprising antigen-binding determinants. Peptide mimetics may also be used. These molecules are usually conformationally restricted organic rings that mimic the structure of a CDR loop and that include antigen-interactive side chains.

20

Synthetic constructs include chimaeric molecules. Thus, for example, humanised (or primatised) antibodies or derivatives thereof are within the scope of the present invention. An example of a humanised antibody is an antibody having human framework regions, but rodent hypervariable regions. Ways of producing chimaeric antibodies are discussed for example by Morrison *et al* in *PNAS*, **81**, 6851-6855 (1984) 25 and by Takeda *et al* in *Nature*. **314**, 452-454 (1985).

Synthetic constructs also include molecules comprising an additional moiety that provides the molecule with some desirable property in addition to antigen binding. For example the moiety may be a label (e.g. a fluorescent or radioactive label). Alternatively, it may be a pharmaceutically active agent.

5

Affibodies are proteins which are found to bind to target proteins with a low dissociation constant. They are selected from phage display libraries expressing a segment of the target protein of interest (Nord K, Gunneriusson E, Ringdahl J, Stahl S, Uhlen M, Nygren PA, Department of Biochemistry and Biotechnology, Royal Institute of Technology (KTH), Stockholm, Sweden).

10

In a further aspect the invention provides an immunogenic composition comprising one or more proteins, polypeptides, peptides, fragments or derivatives thereof, or nucleotide sequences described herein. The immunogenic composition may include nucleic acid sequences ID-65 and/or ID-66 as described herein. Alternatively, the immunogenic composition may comprise proteins/polypeptides including ID-65, ID-83, ID-89, ID-93 and/or ID-96 as described herein, or fragments or derivatives thereof. A composition of this sort may be useful in the treatment or prevention of Group B Streptococcus infection in subject. In a preferred aspect of the invention the immunogenic composition is a vaccine.

15

20

In other aspects the invention provides:

- i) Use of an immunogenic composition as described herein in the preparation of a medicament for the treatment or prophylaxis of Group B Streptococcus infection. Preferably the medicament is a vaccine.

25

- ii) A method of detection of Group B Streptococcus which comprises the step of bringing into contact a sample to be tested with at least one antibody, affibody, or a derivative thereof, as described herein.
- 5      iii) A method of detection of Group B Streptococcus which comprises the step of bringing into contact a sample to be tested with at least one protein, polypeptide, peptide, fragments or derivatives as described herein.
- 10      iv) A method of detection of Group B Streptococcus which comprises the step of bringing into contact a sample to be tested with at least one nucleic acid molecule as described herein.
- v) A kit for the detection of Group B Streptococcus comprising at least one antibody, affibody, or derivatives thereof, described herein.
- 15      vi) A kit for the detection of Group B Streptococcus comprising at least one Group B Streptococcus protein, polypeptide, peptide, fragment or derivative thereof, as described herein.
- 20      vii) A kit for the detection of Group B Streptococcus comprising at least one nucleic acid of the invention.

As described previously, the novel proteins described herein are identified and isolated using a screening method which specifically identifies those Group B  
 25 Streptococcus genes encoding bacterial cell envelope associated or secreted proteins.

Given that the inventors have identified a group of important proteins, such proteins are potential targets for anti-microbial therapy. It is necessary, however, to

determine whether each individual protein is essential for the organism's viability. Thus, the present invention also provides a method of determining whether a protein or polypeptide as described herein represents a potential anti-microbial target which comprises inactivating said protein and determining whether Group B Streptococcus is still viable.

A suitable method for inactivating the protein is to effect selected gene knockouts, ie prevent expression of the protein and determine whether this results in a lethal change. Suitable methods for carrying out such gene knockouts are described in Li *et al* , *P.N.A.S.*, **94**:13251-13256 (1997) and Kolkman *et al.*, *Journal of Biological Chemistry* **272**: 19502-19508 (1997); Kolkman *et al.*, *Journal of Bacteriology* **178**: 3736-3741 (1996).

In a final aspect the present invention provides the use of an agent capable of antagonising, inhibiting or otherwise interfering with the function or expression of a protein or polypeptide of the invention in the manufacture of a medicament for use in the treatment or prophylaxis of Group B Streptococcus infection.

The invention will now be described by means of the following examples which should not in any way be construed as limiting. The examples refer to the figures in which:

Fig 1: (A) Shows a number of full length nucleotide sequences encoding antigenic Group B Streptococcus proteins and the corresponding amino acid sequences.

Fig 2: Shows the results of vaccine trials using the proteins ID-65 and ID-66;

Fig 3: Shows a number of oligonucleotide primers used in the screening process

**nucS1** primer designed to amplify a mature form of the *nuc* A gene

**nucS2-** primer designed to amplify a mature form of the *nuc* A gene.

5 **nucS3** primer designed to amplify a mature form of the *nuc* A gene

**nucR** primer designed to amplify a mature form of the *nuc* A gene

**nucseq** primer designed to sequence DNA cloned into the pTREP-Nuc vector

**pTREPF** nucleic acid sequence containing recognition site for *ECORV*. Used for cloning fragments into pTREX7.

10 **pTREPR** nucleic acid sequence containing recognition site for *BAMH1*. Used for cloning fragments into pTREX7.

**PUCF** forward sequencing primer, enables direct sequencing of cloned DNA fragments.

15 **VR** example of gene specific primer used to obtain further antigen DNA sequence by the method of DNA walking.

**V1** example of gene specific primer used to obtain further antigen DNA sequence by the method of DNA walking.

**V2** example of gene specific primer used to obtain further antigen DNA sequence by the method of DNA walking.

20

Fig 4: (i) Schematic presentation of the nucleotide sequence of the unique gene cloning site immediately upstream of the mature *nuc* gene in pTREP1-*nuc1*, pTREP1-*nuc2* and pTREP1-*nuc3*. Each of the pTREP-*nuc* vectors contain an *EcoRV* (a *SmaI* site in pTREP1-*nuc2*) cleavage site which allows cloning of genomic DNA fragments in 3 different frames with respect to the mature *nuc* gene.

25

(ii) A physical and genetic summary map of the pTREP1-*nuc* vectors. The expression cassette incorporating *nuc*, the macrolides, lincosamides and

streptogramin B (MLS) resistance determinant, and the replicon (rep) *Ori-pAMβ1* are depicted (not drawn to scale).

(iii) Schematic presentation of the expression cassette showing the various sequence elements involved in gene expression and location of unique restriction endonuclease sites (not drawn to scale).

5

Fig 5: SDS-PAGE analysis of a purified preparation of the His-tagged ID-65 and ID-83 protein antigens (predicted molecular weights of 57,144 and 25,000 daltons respectively) on a 12% polyacrylamide gel. Lanes: MW, molecular weight standards; 1, His-tagged ID-65 protein; 2, His-tagged ID-83 protein

10

Fig 6: SDS PAGE analysis of a purified preparation of the His-tagged ID-93 protein antigen (predicted molecular weight = 28,000 daltons) on a 12% polyacrylamide gel.

15

Lanes: MW, molecular weight standards; 1, His-tagged ID-93 protein.

Fig 7: SDS PAGE analysis of a purified preparation of the His-tagged ID-89 and ID-96 protein antigens (predicted molecular weights of 35,000 and 31,000 daltons respectively) on a 12% polyacrylamide gel.

20

Lanes: MW, molecular weight standards; 1, His-tagged ID-89 protein; 2, His-tagged ID-96 protein.

Fig 8: IgG Titres against the ID-65 and ID-83 proteins

25

1 = ID-65 + Alum Group - Bleed at 5 weeks

2 = PBS + Alum Control Group - Bleed at 5 weeks

(For groups 1 and 2, ELISAs were performed on purified ID-65 protein)

3 = ID-83 + Alum Group - Bleed at 5 weeks



17.

4 = PBS + Alum Control Group – Bleed at 5 weeks

(For groups 3 and 4, ELISAs were performed on purified ID-83 protein)

Fig 9: Shows the results of vaccine trials using the protein ID-93.

5

Fig 10: IgG titres against the ID-93 protein.

1 = ID-93 + Alum Group – Bleed at 3 weeks

2 = ID-93 + Alum Group – Bleed at 6 weeks

3 = PBS + Alum Control Group – Bleed at 3 weeks

10

4 = PBS + Alum Control Group – Bleed at 6 weeks

Fig 11: IgG titres against the ID-89 and ID-96 proteins

1 = ID-89 + TitreMax Gold Group – Bleed at 3 weeks

2 = ID-89 + TitreMax Gold – Bleed at 6 weeks

15

3 = PBS + TitreMax Gold Control Group – Bleed at 3 weeks

4 = PBS + TitreMax Gold Control Group – Bleed at 6 weeks

5 = ID-96 + TitreMax Gold Group – Bleed at 3 weeks

6 = ID-96 + TitreMax Gold Group – Bleed at 6 weeks

7 = PBS + TitreMax Gold Control Group – Bleed at 3 weeks

20

8 = PBS + TitreMax Gold Control Group – Bleed at 6 weeks

For Groups 1-4, ELISAs were performed on purified ID-89 protein.

For Groups 5-6, ELISAs were performed on purified ID-96 protein.

Fig 12: Southern blot analysis of genomic DNA. Genomic DNA from each of the strains listed in Table 7 was digested completely with *Hin* DIII (NEB) and electrophoresed at 40 Volts for 6 hours in 0.8% agarose, transferred onto Hybond N<sup>+</sup> (Amersham) membrane by Southern blot and hybridised with the

25

digoxigenin-labelled *rib* gene probe. Specifically bound DNA probe was identified using the DIG Nucleic Acid Detection Kit (Boehringer Mannheim).

5 Fig 13: Southern blot analysis of genomic DNA. Genomic DNA from each of the strains listed in Table 6 was digested completely with *Hin* DIII (NEB) and electrophoresed at 40 Volts for 6 hours in 0.8% agarose, transferred onto Hybond N<sup>+</sup> (Amersham) membrane by Southern blot and hybridised with the digoxigenin-labelled ID-65 gene probe. Specifically bound DNA probe was identified using the DIG Nucleic Acid Detection Kit (Boehringer Mannheim).

10 Fig 14: Southern blot analysis of genomic DNA. Genomic DNA from each of the strains listed in Table 6 was digested completely with *Hin* DIII (NEB) and electrophoresed at 40 Volts for 6 hours in 0.8% agarose, transferred onto Hybond N<sup>+</sup> (Amersham) membrane by Southern blot and hybridised with the digoxigenin-labelled ID-89 gene probe. Specifically bound DNA probe was identified using the DIG Nucleic Acid Detection Kit (Boehringer Mannheim).

15 Fig 15: Southern blot analysis of genomic DNA. Genomic DNA from each of the strains listed in Table 6 was digested completely with *Hin* DIII (NEB) and electrophoresed at 40 Volts for 6 hours in 0.8% agarose, transferred onto Hybond N<sup>+</sup> (Amersham) membrane by Southern blot and hybridised with the digoxigenin-labelled ID-93 gene probe. Specifically bound DNA probe was identified using the DIG Nucleic Acid Detection Kit (Boehringer Mannheim).

20 Fig 16: Southern blot analysis of genomic DNA. Genomic DNA from each of the strains listed in Table 6 was digested completely with *Eco* RI (NEB) and electrophoresed at 40 Volts for 6 hours in 0.8% agarose, transferred onto Hybond N<sup>+</sup> (Amersham) membrane by Southern blot and hybridised with the

digoxigenin-labelled ID-96 gene probe. Specifically bound DNA probe was identified using the DIG Nucleic Acid Detection Kit (Boehringer Mannheim).

5

### Example 1

Gene/partial gene sequences putatively encoding exported proteins in *S. agalactiae* have been identified, unless stated otherwise, using the nuclease screening system described herein *vis*, the LEEP (Lactococcus Expression of Exported Proteins) system. These have been further analysed to remove artefacts. The nucleotide sequences of genes identified using the screening system have been characterised using a number of parameters described below.

1. All putative surface proteins are analysed for leader/signal peptide sequences. Bacterial signal peptide sequences share a common design. They are characterised by a short positively charged N-terminus (N region) immediately preceding a stretch of hydrophobic residues (central portion-h region) followed by a more polar C-terminal portion which contains the cleavage site (c-region). Computer software is used to perform hydropathy profiling of putative proteins (Marcks, *Nuc. Acid. Res.*, **16**:1829-1836 (1988)) which is used to identify the distinctive hydrophobic portion (h-region) typical of leader peptide sequences. In addition, the presence/absence of a potential ribosomal binding site (Shine-Dalgarno sequence required for translation) is also noted.

2. All putative surface protein sequences are used to search the OWL sequence database which includes a translation of the GENBANK and SWISSPROT database.. This allows identification of similar sequences which may have been previously characterised not only at the sequence level but at a functional level. It

may also provide information indicating that these proteins are indeed surface related and not artefacts.

3. Putative *S. agalactiae* surface proteins are also assessed for their novelty. Some of the identified proteins may or may not possess a typical leader peptide sequence and may not show homology with any DNA/protein sequences in the database. Indeed these proteins may indicate the primary advantage of our screening method, i.e. isolating atypical surface-related proteins, which would have been missed in all previously described screening protocols.

The construction of three reporter vectors and their use in *L. lactis* to identify and isolate genomic DNA fragments from pathogenic bacteria encoding secreted or surface associated proteins is now described.

#### Construction of the pTREP1-*nuc* series of reporter vectors

##### (a) Construction of expression plasmid pTREP1

The pTREP1 plasmid is a high-copy number (40-80 per cell) theta-replicating gram positive plasmid, which is a derivative of the pTREX plasmid which is itself a derivative of the previously published pIL253 plasmid. pIL253 incorporates the broad Gram-positive host range replicon of pAM $\beta$ 1 (Simon and Chopin, *Biochemie* 70: 559-566 (1988)) *L. lactis* sex-factor. pIL253 also lacks the *tra* function which is necessary for transfer or efficient mobilisation by conjugative parent plasmids exemplified by pIL501. The Enterococcal pAM $\beta$ 1 replicon has previously been transferred to various species including *Streptococcus*, *Lactobacillus* and *Bacillus* species as well as *Clostridium acetobutylicum*, (LeBlanc *et al.*, *Proceedings of the National Academy of Science USA* 75:3484-3487 (1978)) indicating the potential broad host range utility. The pTREP1 plasmid represents a constitutive transcription vector.

The pTREX vector was constructed as follows. An artificial DNA fragment containing a putative RNA stabilising sequence, a translation initiation region (TIR), a multiple cloning site for insertion of the target genes and a transcription terminator was created by annealing 2 complementary oligonucleotides and extending with Tfl DNA polymerase. The sense and anti-sense oligonucleotides contained the recognition sites for NheI and BamHI at their 5' ends respectively to facilitate cloning. This fragment was cloned between the XbaI and BamHI sites in pUC19NT7, a derivative of pUC19 which contains the T7 expression cassette from pLET1 (Wells *et al.*, *J. Appl. Bacteriol.* 74:629-636 (1993)) cloned between the EcoRI and HindIII sites. The resulting construct was designated pUCLEX. The complete expression cassette of pUCLEX was then removed by cutting with HindIII and blunting followed by cutting with EcoRI before cloning into EcoRI and SacI (blunted) sites of pIL253 to generate the vector pTREX (Wells and Schofield, *In Current advances in metabolism, genetics and applications-NATO ASI Series. H* 98:37-62. (1996)). The putative RNA stabilising sequence and TIR are derived from the *Escherichia coli* T7 bacteriophage sequence and modified at one nucleotide position to enhance the complementarity of the Shine Dalgarno (SD) motif to the ribosomal 16s RNA of *Lactococcus lactis* (Schofield *et al.* pers. coms. University of Cambridge Dept. Pathology.).

A *Lactococcus lactis* MG1363 chromosomal DNA fragment exhibiting promoter activity which was subsequently designated P7 was cloned between the EcoRI and BglII sites present in the expression cassette, creating pTREX7. This active promoter region had been previously isolated using the promoter probe vector pSB292 (Waterfield *et al.*, *Gene* 165:9-15 (1995)). The promoter fragment was amplified by PCR using the Vent DNA polymerase according to the manufacturer.

The pTREP1 vector was then constructed as follows. An artificial DNA fragment which included a transcription terminator, the forward pUC sequencing primer, a promoter multiple cloning site region and a universal translation stop sequence was created by annealing two overlapping partially complementary synthetic  
5 oligonucleotides together and extending with sequenase according to manufacturers instructions. The sense and anti-sense (pTREP<sub>F</sub> and pTREP<sub>R</sub>) oligonucleotides contained the recognition sites for EcoRV and BamHI at their 5' ends respectively to facilitate cloning into pTREX7. The transcription terminator was that of the *Bacillus penicillinase* gene, which has been shown to be effective in *Lactococcus* (Jos *et al.*,  
10 *Applied and Environmental Microbiology* 50:540-542 (1985)). This was considered necessary as expression of target genes in the pTREX vectors was observed to be leaky and is thought to be the result of cryptic promoter activity in the origin region (Schofield *et al.* pers. coms. University of Cambridge Dept. Pathology.). The forward pUC primer sequencing was included to enable direct sequencing of cloned  
15 DNA fragments. The translation stop sequence which encodes a stop codon in 3 different frames was included to prevent translational fusions between vector genes and cloned DNA fragments. The pTREX7 vector was first digested with EcoRI and blunted using the 5' - 3' polymerase activity of T4 DNA polymerase (NEB) according to manufacturer's instructions. The EcoRI digested and blunt ended  
20 pTREX7 vector was then digested with Bgl II thus removing the P7 promoter. The artificial DNA fragment derived from the annealed synthetic oligonucleotides was then digested with EcoRV and Bam HI and cloned into the EcoRI(blunted)-Bgl II digested pTREX7 vector to generate pTREP. A *Lactococcus lactis* MG1363 chromosomal promoter designated P1 was then cloned between the EcoRI and BglII  
25 sites present in the pTREP expression cassette forming pTREP1. This promoter was also isolated using the promoter probe vector pSB292 and characterised by Waterfield *et al.*, (1995) [*supra*]. The P1 promoter fragment was originally amplified by PCR using vent DNA polymerase according to manufacturers

instructions and cloned into the pTREX as an EcoRI-BglII DNA fragment. The EcoRI-BglII P1 promoter containing fragment was removed from pTREX1 by restriction enzyme digestion and used for cloning into pTREP (Schofield *et al.* pers. coms. University of Cambridge, Dept. Pathology.).

5

**(b) PCR amplification of the *S. aureus* nuc gene.**

The nucleotide sequence of the *S. aureus* nuc gene (EMBL database accession number V01281) was used to design synthetic oligonucleotide primers for PCR  
10 amplification. The primers were designed to amplify the mature form of the nuc gene designated nucA which is generated by proteolytic cleavage of the N-terminal 19 to 21 amino acids of the secreted propeptide designated Snase B (Shortle, 1983 [supra]). Three sense primers (nucS1, nucS2 and nucS3, shown in figure 3) were designed, each one having a blunt-ended restriction endonuclease cleavage site for  
15 EcoRV or SmaI in a different reading frame with respect to the nuc gene. Additionally BglII and BamHI were incorporated at the 5' ends of the sense and anti-sense primers respectively to facilitate cloning into BamHI and BglII cut pTREP1. The sequences of all the primers are given in figure 3. Three nuc gene DNA fragments encoding the mature form of the nuclease gene (NucA) were amplified by  
20 PCR using each of the sense primers combined with the anti-sense primer. The nuc gene fragments were amplified by PCR using *S. aureus* genomic DNA template, Vent DNA Polymerase (NEB) and the conditions recommended by the manufacturer. An initial denaturation step at 93°C for 2 min was followed by 30 cycles of denaturation at 93°C for 45 sec, annealing at 50°C for 45 seconds, and extension at  
25 73°C for 1 minute and then a final 5 min extension step at 73°C. The PCR amplified products were purified using a Wizard clean up column (Promega) to remove unincorporated nucleotides and primers.

**(c) Construction of the pTREP1-*nuc* vectors**

The purified *nuc* gene fragments described in section b were digested with Bgl II and BamHI using standard conditions and ligated to BamHI and BglII cut and dephosphorylated pTREP1 to generate the pTREP1-*nuc*1, pTREP1-*nuc*2 and pTREP1-*nuc*3 series of reporter vectors. These vectors are described in figure 4. General molecular biology techniques were carried out using the reagents and buffers supplied by the manufacturer or using standard techniques (Sambrook and Maniatis, Molecular cloning: A laboratory manual. Cold Spring Harbor Laboratory Press: Cold Spring Harbour (1989)). In each of the pTREP1-*nuc* vectors the expression cassette comprises a transcription terminator, lactococcal promoter P1, unique cloning sites (Bgl II, EcoRV or SmaI) followed by the mature form of the *nuc* gene and a second transcription terminator. Note that the sequences required for translation and secretion of the *nuc* gene were deliberately excluded in this construction. Such elements can only be provided by appropriately digested foreign DNA fragments (representing the target bacterium) which can be cloned into the unique restriction sites present immediately upstream of the *nuc* gene.

**(d) Screening for secreted proteins in Group B Streptococcus.**

Genomic DNA isolated from Group B Streptococcus (*S. agalactiae*) was digested with the restriction enzyme Tru9I. This enzyme which recognises the sequence 5'-TTAA -3' was used because it cuts A/T rich genomes efficiently and can generate random genomic DNA fragments within the preferred size range (usually averaging 0.5 - 1.0 kb). This size range was preferred because there is an increased probability that the P1 promoter can be utilised to transcribe a novel gene sequence. However, the P1 promoter may not be necessary in all cases as it is possible that many Streptococcal promoters are recognised in *L. lactis*. DNA fragments of different size ranges were purified from partial Tru9I digests of *S. agalactiae* genomic DNA. As



the Tru 9I restriction enzyme generates staggered ends the DNA fragments had to be made blunt ended before ligation to the EcoRV or SmaI cut pTREP1-*nuc* vectors. This was achieved by the partial fill-in enzyme reaction using the 5'-3' polymerase activity of Klenow enzyme. Briefly Tru9I digested DNA was dissolved in a solution (usually between 10-20  $\mu$ l in total) supplemented with T4 DNA ligase buffer (New England Biolabs; NEB) (1X) and 33  $\mu$ M of each of the required dNTPs, in this case dATP and dTTP. Klenow enzyme was added (1 unit Klenow enzyme (NEB) per  $\mu$ g of DNA) and the reaction incubated at 25°C for 15 minutes. The reaction was stopped by incubating the mix at 75°C for 20 minutes. EcoRV or SmaI digested pTREP-*nuc* plasmid DNA was then added (usually between 200-400 ng). The mix was then supplemented with 400 units of T4 DNA ligase (NEB) and T4 DNA ligase buffer (1X) and incubated overnight at 16°C. The ligation mix was precipitated directly in 100% Ethanol and 1/10 volume of 3M sodium acetate (pH 5.2) and used to transform *L. lactis* MG1363 (Gasson, *J. Bacteriol.* 154:1-9 (1983)). Alternatively, the gene cloning site of the pTREP-*nuc* vectors also contains a BglII site which can be used to clone for example Sau3AI digested genomic DNA fragments.

*L. lactis* transformant colonies were grown on brain heart infusion agar and nuclease secreting (*Nuc*<sup>+</sup>) clones were detected by a toluidine blue-DNA-agar overlay (0.05 M Tris pH 9.0, 10 g of agar per litre, 10 g of NaCl per liter, 0.1 mM CaCl<sub>2</sub>, 0.03 % wt/vol. salmon sperm DNA and 90 mg of Toluidine blue O dye) essentially as described by Shortle, 1983 [*supra*], and Le Loir *et al.*, 1994 [*supra*]). The plates were then incubated at 37°C for up to 2 hours. Nuclease secreting clones develop an easily identifiable pink halo. Plasmid DNA was isolated from *Nuc*<sup>+</sup> recombinant *L. lactis* clones and DNA inserts were sequenced on one strand using the *NucSeq* sequencing primer described in figure 3, which sequences directly through the DNA insert.

## Example 2

### Preparation of a *S. agalactiae* standard inoculum

#### 5 Strain validation

*S. agalactiae* serotype III (strain 97/0099) is a recent clinical isolate derived from the cerebral spinal fluid of a new born baby suffering from meningitis. This haemolytic strain of Group B Streptococcus was epidemiologically tested and validated at the Respiratory and Systemic Infection Laboratory, PHLS Central Public Health Laboratory, 61 Colindale Avenue, London NW9 5HT. The strain was subcultured only twice prior to its arrival in the laboratory. Upon its arrival on an agar slope, a sweep of 4-5 colonies was immediately used to inoculate a Todd Hewitt/5% horse blood broth which was incubated overnight statically at 37°C. 0.5 ml aliquots of this overnight culture were then used to make 20% glycerol stocks of the bacterium for long-term storage at -70°C. Glycerol stocks were streaked on Todd Hewitt/5% horse blood agar plates to confirm viability.

#### *In vivo* passaging of Group B Streptococcus

A frozen culture (described under strain validation) of *S. agalactiae* serotype III (strain 97/0099) was streaked to single colonies on Todd-Hewitt/5% blood agar plates, which were incubated overnight at 37°C. A sweep of 4-5 colonies was used to inoculate a Todd Hewitt/5% horse blood broth, which was again incubated overnight. A 0.5 ml aliquot from this overnight culture was used to inoculate a 50 ml Todd Hewitt broth (1:100 dilution) which was incubated at 37°C. 10-fold serial dilutions of the overnight culture were made (since virulence of this strain was unknown) and each was passaged intra-peritoneally (IP) in CBA/ca mice in duplicate. Viable counts were performed on the various inocula used in the passage. Groups of mice were challenged with various concentrations of the pathogen ranging from  $10^8$  to  $10^4$  colony forming units (cfu). Mice that developed symptoms were terminally anaesthetized and cardiac punctures were performed (Only mice that had

been challenged with the highest doses, i.e.  $1 \times 10^8$  cfu, developed symptoms). The retrieved unclotted blood was used to inoculate directly a 50ml serum broth (Todd Hewitt/20% inactivated foetal calf serum). The culture was constantly monitored and allowed to grow to late logarithmic phase. The presence of blood in the medium  
5 interfered with OD<sub>600nm</sub> readings as it was being increasingly lysed with increasing growth of the bacterium, hence the requirement to constantly monitor the culture. Upon reaching late logarithmic phase/early stationary phase, the culture was transferred to a fresh 50 ml tube in order to exclude dead bacterial cells and remaining blood cells which would have sedimented at the bottom of the tube. 0.5  
10 ml aliquots were then transferred to sterile cryovials, frozen in liquid nitrogen and stored at  $-70^{\circ}\text{C}$ . A viable count was carried out on a single standard inoculum aliquot in order to determine bacterial numbers. This was determined to be approximately  $5 \times 10^8$  cfu per ml.

15 **Intra-peritoneal Challenge and virulence testing of Group B Streptococcus standard inoculum**

To determine if the standard inoculum was suitably virulent for use in a vaccine trial, challenges were carried out using a dose range. Frozen standard inoculum strain aliquots were allowed to thaw at room temperature. From viable count data the  
20 number of cfu per ml was already known for the standard inoculum. Initially, serial dilutions of the standard inoculum were made in Todd Hewitt broth and mice were challenged intra-peritoneally with doses ranging from  $1 \times 10^8$  to  $1 \times 10^4$  cfu in a  $500 \mu\text{l}$  volume of Todd Hewitt broth. The survival times of mouse groups injected with different doses of the bacterium were compared. The standard inoculum was  
25 determined to be suitably virulent and a dose of  $1 \times 10^6$  cfu was considered close to optimal for further use in vaccine trials. Further optimisation was carried out by comparing mice challenged with doses ranging between  $5 \times 10^5$  and  $5 \times 10^6$  cfu. The optimal dose was estimated to be approximately  $2.5 \times 10^6$  cfu. This represented

a 100% lethal dose and was repeatedly consistent with end-points as determined by survival times being clustered within a narrow time-range. Throughout all these experiments, challenged mice were constantly monitored to clarify symptoms, stages of symptom development as well as calculating survival times.

5

**Screening Group B Streptococcal LEEP derived genes in DNA vaccination experiments.**

**pcDNA3.1+ as a DNA vaccine vector**

10 The commercially available pcDNA3.1+ plasmid (Invitrogen), referred to as pcDNA3.1 henceforth, was used as a vector in all DNA immunisation experiments involving gene targets derived using the LEEP system unless stated otherwise. pcDNA 3.1 is designed for high-level stable and transient expression in mammalian cells and has been used widely and successfully as a host vector to test candidate  
15 genes from a variety of pathogens in DNA vaccination experiments (Zhang *et al.*, *Infection and Immunity* 176: 1035-40 (1997); Kurar and Splitter, *Vaccine* 15: 1851-57 (1997); Anderson *et al.*, *Infection and Immunity* 64: 3168-3173 (1996)).

20 The vector possesses a multiple cloning site which facilitates the cloning of multiple gene targets downstream of the human cytomegalovirus (CMV) immediate-early promoter/enhancer which permits efficient, high-level expression of the target gene in a wide variety of mammalian cells and cell types including both muscle and immune cells. This is important for optimal immune response as it remains unknown  
25 as to which cells types are most important in generating a protective response *in vivo*. The plasmid also contains the ColE1 origin of replication which allows convenient high-copy number replication and growth in *E. coli* and the ampicillin resistance gene (B- lactamase) for selection in *E. coli*. In addition pcDNA 3.1

possesses a T7 promoter/priming site upstream of the MCS which allows for *in vitro* transcription of a cloned gene in the sense orientation.

#### **Preparation of DNA vaccines**

5 Oligonucleotide primers were designed for each individual gene of interest derived using the LEEP system unless stated otherwise. Each gene was examined thoroughly, and where possible, primers were designed such that they targeted that portion of the gene believed to encode only the mature portion of the protein (APPENDIX I); the intention being to express those sequences that encode only the  
10 mature portion of a target gene protein to would facilitate its correct folding when expressed in mammalian cells. For example, in the majority of cases primers were designed such that putative N-terminal signal peptide sequences would not be included in the final amplification product to be cloned into the pcDNA3.1 expression vector. The signal peptide directs the polypeptide precursor to the cell  
15 membrane via the protein export pathway where it is normally cleaved off by signal peptidase I (or signal peptidase II if a lipoprotein). Hence the signal peptide does not make up any part of the mature protein whether it be displayed on the bacterium's surface or secreted. Where an N-terminal leader peptide sequence was not immediately obvious, primers were designed to target the whole of the gene  
20 sequence for cloning and ultimately, expression in pcDNA3.1.

All forward and reverse oligonucleotide primers incorporated appropriate restriction enzyme sites to facilitate cloning into the pcDNA3.1 MCS region. All forward  
25 primers were also designed to include the conserved Kozak nucleotide sequence 5'-gccacc-3' immediately upstream of an 'atg' translation initiation codon in frame with the target gene insert. The Kozak sequence facilitates the recognition of initiator sequences by eukaryotic ribosomes. Typically, a forward primer incorporating a BamH1 restriction enzyme site the primer would begin with the sequence 5'-

cgggatccgccaccatg-3', followed by a sequence homologous to the 5' end of that part of a gene being amplified. All reverse primers incorporated a Not I restriction enzyme site sequence 5' -ttcggccgc-3'. All gene-specific forward and reverse primers were designed with compatible melting temperatures to facilitate their amplification.

All gene targets were amplified by PCR from *S. agalactiae* genomic DNA template using Vent DNA polymerase (NEB) or rTth DNA polymerase (PE Applied Biosystems) using conditions recommended by the manufacturer. A typical amplification reaction involved an initial denaturation step at 95°C for 2 minutes followed by 35 cycles of denaturation at 95°C for 30 seconds, annealing at the appropriate melting temperature for 30 seconds, and extension at 72°C for 1 minute (1 minute per kilobase of DNA being amplified). This was followed by a final extension period at 72°C for 10 minutes. All PCR amplified products were extracted once with phenol chloroform (2:1:1) and once with chloroform (1:1) and ethanol precipitated. Specific DNA fragments were isolated from agarose gels using the QIAquick Gel Extraction Kit (Qiagen). The purified amplification gene DNA fragments were digested with the appropriate restriction enzymes and cloned into the pcDNA3.1 plasmid vector using *E. coli* as a host. Successful cloning and maintenance of genes was confirmed by restriction mapping and by DNA sequencing. Recombinant plasmid DNA was isolated on a large scale (> 1.5 mg) using Plasmid Mega Kits (Qiagen).

#### DNA vaccination trials

DNA vaccine trials in mice were accomplished by the administration of DNA to 6 week old CBA/ca mice (Harlan, UK). Mice to be vaccinated were divided into groups of six and each group was immunised with recombinant pcDNA3.1 plasmid DNA containing a specific target-gene sequence derived using the LEEP system unless stated otherwise. A total of 100 µg of DNA in Dulbecco's PBS (Sigma) was

injected intramuscularly into the tibialis anterior muscle of both hind legs. Four weeks later this procedure was repeated using the same amount of DNA. For comparison, control mice groups were included in all vaccine trials. These control groups were either not DNA-vaccinated or were immunised with non-recombinant pcDNA3.1 plasmid DNA only, using the same time course described above. Four weeks after the second immunisation, all mice groups were challenged intra-peritoneally with a lethal dose of *S. agalactiae* serotype III (strain 97/0099). The actual number of bacteria administered was determined by plating serial dilutions of the inoculum on Todd-Hewitt/5% blood agar plates. All mice were killed 3 or 4 days after infection. During the infection process, challenged mice were monitored for the development of symptoms associated with the onset of *S. agalactiae* induced-disease. Typical symptoms in an appropriate order included piloerection, an increasingly hunched posture, discharge from eyes, increased lethargy and reluctance to move which was often the result of apparent paralysis in the lower body/hind leg region. The latter symptoms usually coincided with the development of a moribund state at which stage the mice were culled to prevent further suffering. These mice were deemed to be very close to death, and the time of culling was used to determine a survival time for statistical analysis. Where mice were found dead, a survival time was calculated by averaging the time when a particular mouse was last observed alive and the time when found dead, in order to determine a more accurate time of death. The results of this trial are shown in Table 1 and presented graphically in Figure 2.

### Interpretation of Results

A positive result was taken as any DNA sequence that was cloned and used in challenge experiments as described above and gave protection against that challenge. DNA sequences were determined to be protective;

-if that DNA sequence gave statistically significant protection to mice as compared to control mice (to a 95% confidence level ( $p > 0.05$ ) as determined using the Mann-Whitney U test .

- 5 -if that DNA sequence was marginal or non-significant using Mann-Whitney but showed some protective features. For example, one or more outlying mice may survive for significantly longer time periods when compared with control mice. Alternatively, the time to first death may also be prolonged when compared to counterpart mice in control groups. It is acceptable to allow marginal or non-significant results to be considered as potential positives when it is possible that the
- 10 clarity of some results may be affected by problems associated with the administration of the DNA vaccine. Indeed, much varied survival times may reflect different levels of immune response between different members of a given group.

Table 1

15 **LEEP DNA immunisation and GBS challenge Experiment**

**Statistical analysis of survival times**

	Mean Survival Times (hours)		
	UnVacc	3-60(ID-65)	3-5(ID-66)
1	27.583	54.416	42.916
2	27.583	31.000	42.916
3	24.583	43.000	32.874
4	22.250	34.916	42.916
5	35.916	38.958	27.333
6	22.250	34.916	30.916
<b>Mean</b>	<b>27.583</b>	<b>40.458</b>	<b>37.791</b>
<b>sd</b>	<b>5.1691</b>	<b>8.9959</b>	<b>7.2860</b>
<b>p value</b>		<b>0.0098</b>	<b>0.0215</b>

- 20 **p value** refers to statistical significance when compared to unvaccinated controls.



**Comment**

## ID-65 (3-60)

Mice immunised with the '3-60 (ID-65)' DNA vaccine exhibited significantly longer survival times when compared with the unvaccinated control group.

**ID-66 (3-5)**

Mice immunised with the '3-5 (ID-66)' DNA vaccine exhibited significantly longer survival times when compared with the unvaccinated control group.

### Example 3

### Expression and Screening Group B Streptococcal LEEP derived Proteins in Protein vaccination experiments.

### Expression of proteins

Prioritised genes ie, those selected on the basis of predicted expression features as deduced from sequence characteristics (as described in Figure 1), were cloned and expressed as recombinant proteins using the pET system (Novagen, Inc., Madison, WI) utilising *Escherichia coli* as a host. Target genes were cloned into the pET28b(+) plasmid expression vector. The pET28b(+) vector is designed for high level expression and purification of target proteins. This vector carries a T7 promoter for transcription of a target gene, followed by an N-terminal His•Tag®/thrombin/T7•Tag® configuration, a multi-cloning site containing unique restriction enzyme sites for cloning purposes, and an optional C-terminal His•Tag sequence. The vector also carries a kanamycin resistance gene for selection purposes and for maintaining target gene expression (pET System Manual, 8<sup>th</sup> edition, Novagen).

### Preparation of protein vaccines

Oligonucleotide primers were designed for each individual target gene derived using the LEEP system unless stated otherwise. Each gene was examined thoroughly. Where possible primers were designed so that they would target that part of the gene predicted to encode only the mature portion of the protein (APPENDIX II). It is hoped that expressing those corresponding to the predicted mature protein only, might facilitate its correct folding when finally expressed *in vitro*. Oligonucleotide primers were designed so that sequences, encoding the putative N-terminal signal peptide of the target protein, would not be included in the final amplification product to be cloned pET28b(+). The signal peptide directs the polypeptide precursor to the cell membrane via the protein export pathway where it is normally cleaved off by signal peptidase I (or signal peptidase II if a lipoprotein). Hence the signal peptide would not be expected to form any part of the mature target protein, whether it be displayed on the bacterium's surface or secreted. For this purpose, classical signal peptides and their cleavage sites were predicted using the DNA Strider™ Program (CEA, France) and the SignalP V1.1 program, which predicts the presence and location of signal peptide cleavage sites in amino acid sequences from different organisms (Nielsen *et al.*, *Protein Engineering* 10: 1-6 (1997)). Where a N-terminal leader peptide sequence was not obvious, primers were designed to include the whole of the gene sequence for cloning and expression.

All oligonucleotide primers were designed to incorporate appropriate restriction enzyme sites to facilitate cloning into the pcDNA3.1 MCS region (APPENDIX II). Forward primers included an *Nco* I (5'-ccatgg-3') or *Nhe* I (5'-gctagc-3') restriction enzyme site and an 'ATG' start codon in-frame with the target gene open reading frame (orf). All reverse primers included a *Not* I restriction enzyme site 5' - gcggtccgc-3' and were designed so that the target gene could be expressed in frame with the C-terminal His•Tag (i.e. the stop codon of the target gene was not

included). Using the *Nco* I and *Not* I, allowed the removal of the N-terminal His•Tag®, thrombin and T7•Tag® DNA sequences. At the same time target genes were cloned immediately downstream of a highly efficient ribosome binding site (from the phage T7 major capsid protein), to facilitate high level expression/translation of the target gene by T7 RNA polymerase, and subsequent purification by means of the C-terminal His•Tag. All target gene-specific forward and reverse primers were designed with compatible melting temperatures to facilitate their amplification.

All gene targets were amplified by PCR from *S. agalactiae* genomic DNA template using Vent DNA polymerase (NEB) using conditions recommended by the manufacturer. A typical amplification reaction involved an initial denaturation step at 95°C for 2 minutes followed by 35 cycles of denaturation at 95°C for 30 seconds, annealing at the appropriate melting temperature for 30 seconds, and extension at 72°C for 1 minute (1 minute per kilobase of DNA being amplified). This was followed by a final extension period at 72°C for 10 minutes. All PCR amplified products were extracted once with phenol:chloroform (2:1:1) and once with chloroform (1:1) and ethanol precipitated. Specific DNA fragments were isolated from agarose gels using the QIAquick Gel Extraction Kit (Qiagen). Purified target gene DNA amplicons were then digested *Nco* I (or *Nhe* I) and *Not* I restriction enzymes, and cloned into *Nco* I and *Not* I digested pET28b(+) plasmid vector using *E. coli* DH5α or *E. coli* BL21 (DE3) as a host. Successful cloning and maintenance of genes was confirmed by restriction mapping.

#### **Determination of target protein expression and solubility**

Glycerol stocks of *E. coli* BL21 DE3 pET28b(+) strains expressing recombinant proteins were used to inoculate 10 ml Luria broth containing Kanamycin (30 µg/ml ) which were grown overnight at 37°C with vigorous shaking (300 rpm).

A 20-40 ml Luria broth containing Kanamycin (30  $\mu$ g/ml) was inoculated with 1:100 dilution of the overnight culture from step 1 and grown at 37°C with vigorous shaking (300 rpm). When the culture reached an OD<sub>600</sub> of between 0.6 and 1.0, IPTG was added to a final concentration of 1mM. Typically cultures were induced for 3 hours. Cells were then harvested by centrifugation at 7000 g for 10 min. The cell pellet was then resuspended in 1/10 volume of lysis buffer (50mM NaH<sub>2</sub>PO<sub>4</sub>, pH.8.0; 300mM NaCl; 10mM imidazole; 10% glycerol). Lysozyme was then added to a final concentration of 1mg/ml, and the suspension was incubated on ice for 30 min. The suspension was then sonicated on ice (six 10-sec bursts at 200-300 W with a 10-sec cooling period. The lysate was then centrifuged at 10,000g for 20 min. The supernatant (containing soluble protein) was transferred to a sterile 2 ml eppendorf. The pellet was resuspended in 2 ml of solubilisation buffer (8 M Urea; 50mM NaH<sub>2</sub>PO<sub>4</sub>, pH.8.0; 300mM NaCl; 10% glycerol). This suspension contained the insoluble protein fraction. Aliquots from both the soluble and insoluble fractions were transferred to new eppendorfs. The protein samples were denatured by adding an equal volume of 2x SDS-PAGE buffer and heating at 95°C for 5 min. Denatured extract samples were then analysed by SDS-PAGE to determine target gene expression and solubility.

### Large scale expression of recombinant target proteins

Glycerol stocks of *E. coli* BL21 DE3 pet28b(+) strains expressing recombinant proteins were used to inoculate 10 ml Luria broth containing Kanamycin ( 30  $\mu$ g/ml ) which were grown overnight at 37°C with vigorous shaking (300 rpm). 5 ml of an overnight culture of a recombinant strain was used to inoculate a 250 ml Luria broth containing kanamycin (30  $\mu$ g/ml) which was grown at 37°C with vigorous shaking (300 rpm). When the culture reached an OD<sub>600</sub> of between 0.6 and 1.0, IPTG was added to a final concentration of 1mM. Typically, cultures were induced for 3 hours. Cultures were then centrifuged to a pellet and stored frozen at -20°C.

### Purification of target antigens.

Ni-NTA agarose (Qiagen LTD, West Sussex, UK; Cat. No. 30210) was used to purify the His-Tagged recombinant proteins. The 6xHis affinity tag which was expressed in frame with the target proteins in pET28b(+), facilitates binding to Ni-NTA. Ni-NTA offers high binding capacity (with minimal non-specific binding) and can bind 5-10 mg of 6xHis-tagged protein per ml of resin. The 6xHis-tag is poorly immunogenic, and at pH 8.0, the tag is small, uncharged and therefore does not generally interfere with the structure and function of the protein (The *QIAexpressionist*, Qiagen Handbook, March 1999).

NOTE: All the proteins (LEEP-derived, unless stated otherwise) described here were purified under denaturing conditions except ID-65. ID-65 was prepared and purified under native conditions.

### Purification under native conditions

The frozen pellet was allowed to thaw on ice for 15 minutes and then resuspended in 10 ml of lysis buffer (50mM NaH<sub>2</sub>PO<sub>4</sub>, pH.8.0; 300mM NaCl;10mM imidazole;

10% glycerol). Lysozyme was then added to a final concentration of 1mg/ml, and the suspension was incubated on ice for 30 min. The suspension was then sonicated on ice (six 10-sec bursts at 200-300 W with a 10-sec cooling period). Dnase I (5  $\mu$ g/ml) was then added to the lysate, which was then incubated on ice for 10-15 min.

5 The lysate was then centrifuged at 10,000 rpm for 20 min at 4°C to pellet cell debris. The clear lysate supernatant was then loaded into a polypropylene column (Qiagen; Cat. No. 34964), bottom cap attached. 1.5 ml of 50% Ni-NTA was then added, the column sealed and the suspension was allowed to mix gently using a rotating wheel for 1-2 hours at 4°C. The column containing the lysate/Ni-NTA mix was then

10 placed upright using a retort stand, and the Ni-NTA was allowed to settle. The bottom cap was removed and the lysate was allowed to flow through. The column was then washed with three to six 4 ml volumes of wash buffer (50mM NaH<sub>2</sub>PO<sub>4</sub>, pH.8.0; 300mM NaCl; 20mM imidazole; 10% glycerol). The protein was then eluted in 0.5 ml aliquots of elution buffer (50mM NaH<sub>2</sub>PO<sub>4</sub>, pH.8.0; 300mM

15 NaCl; 500mM imidazole; 10% glycerol). Eluate fractions were then analysed by SDS-PAGE and those containing the protein were pooled and dialysed against a PBS (pH 7.0)-glycerol (10%) solution.

#### **Purification and refolding under denaturing conditions**

20 The frozen pellet was allowed to thaw on ice for 15 minutes and then resuspended in 10 ml of buffer containing 8 M Urea, 300 mM NaCl, 10% glycerol, 0.1 M NaH<sub>2</sub>PO<sub>4</sub>, pH.8.0, and 10 mM imidazole. The cells were then lysed by gentle vortexing for 1 hour at room temperature. The lysate was then centrifuged at

25 10,000g for 20 minutes to pellet cellular debris. The clear lysate supernatant was then loaded into a polypropylene column (Qiagen; Cat. No. 34964), bottom cap attached. 1.5 ml of 50% Ni-NTA slurry was then added, the column sealed and the suspension was allowed to mix gently using a rotating wheel for 1-2 hours at room

temperature. The column containing the lysate/Ni-NTA mix was then placed upright using a retort stand, and the Ni-NTA was allowed to settle. The bottom cap was removed and the lysate was allowed to flow through. The column was then washed with 4-8 ml of buffer containing 8 M Urea, 300 mM NaCl, 10% glycerol, 0.1 M NaH<sub>2</sub>PO<sub>4</sub>, pH 8.0, and 10 mM imidazole. The resin was then washed with a gradient of 6 to 0 M in a buffer containing 0.1 M NaH<sub>2</sub>PO<sub>4</sub>, pH 8.0, 300 mM NaCl and 10% glycerol to facilitate the slow removal of urea and gradual refolding of target protein. The resin was then washed with a buffer containing 0.1 M NaH<sub>2</sub>PO<sub>4</sub>, pH 7.0, 500 mM NaCl and 10% glycerol. The recombinant protein was then eluted in 0.5 ml aliquots with 500 mM Imidazole in 0.1 M NaH<sub>2</sub>PO<sub>4</sub>, pH 7.0, 500 mM NaCl and 10% glycerol. The fractions were analysed on SDS-PAGE and those containing the protein were pooled and dialysed against a PBS (pH 7.0)-glycerol (10%) solution.

All purified proteins were analysed by SDS-PAGE, as shown in Figures 5, 6 and 7, prior to their use as antigens in immunisation and vaccination experiments.

### Protein Vaccinations

Vaccines were composed of the target protein in phosphate buffered saline/10% glycerol and mixed with aluminium hydroxide (alum) (Imject<sup>®</sup> Alum, Pierce, Rockford, Ill.). Each dose (unless otherwise stated) of vaccine contained 25 µg of purified protein in 50 µl of PBS/10% glycerol, mixed with 50 µl of alum. Groups of 6-8 CBA/ca mice (Harlan, UK) were immunised subcutaneously with the vaccines and again 4 weeks later. A control group received 100 µl dose of PBS/10% glycerol with alum. All vaccinated groups consisted of 6 mice. Mice were challenged at 7 weeks (unless otherwise stated). Mice were injected intraperitoneally (i.p.) with between 2.5-5 X 10<sup>6</sup> bacteria diluted in 0.5 ml Todd-Hewitt broth. Deaths were recorded daily for 7 days. The challenged mice were observed daily for signs of illness. Typical symptoms in an appropriate order included piloerection, an

increasingly hunched posture, discharge from eyes, increased lethargy and reluctance to move which was often the result of apparent paralysis in the lower body/hind leg region. The latter symptoms usually coincided with the development of a moribund state at which stage the mice were culled to prevent further suffering. These mice were deemed to be very close to death, and the time of culling was used to determine a survival time for statistical analysis. Where mice were found dead, a survival time was calculated by averaging the time when a particular mouse was last observed alive and the time when found dead, in order to determine a more accurate time of death.

#### **Analysis of antibody responses**

Mice (6 per group) were immunised with two doses of vaccine with a four week interval. Mice were tail bled at 3 weeks and 6 weeks post primary vaccination to obtain sera. Total Immunoglobulin G (IgG) titres to the vaccine protein component in sera were determined by enzyme-linked immunosorbent assay (ELISA), using the original purified protein as the coating antigen.

#### **Standard ELISA protocol**

##### **Solutions**

##### **Carbonate/bicarbonate buffer, pH 9.8**

0.80g  $\text{Na}_2\text{CO}_3$

1.46g  $\text{NaHCO}_3$

pH to 9.6 using HCl

Add distilled water ( $\text{dH}_2\text{O}$ ) to a final volume of 500ml.

#### **p-NITROPHENYL PHOSPHATE SUBSTRATE**



# **Diethanolamine Buffer, pH 9.8**

48.5 ml diethanolamine

pH to 9.8 using 1M HCl

Add dH<sub>2</sub>O to a final volume of 500ml

5

NOTE: ELISAs were optimised for each protein submitted for immunisation.

## **PROTOCOL**

- 10 1. ELISA plates (Greiner labortechnik 96 well plates: Cat. No. 655061) with an appropriate concentration of recombinant protein diluted in carbonate/bicarbonate buffer (50 µl/well). Cover plates with plastic or foil and leave overnight at 4°C.
2. Quickly wash plates twice in a tub/container containing PBS/0.05% Tween-20 and then pat dry.
- 15 3. Block plates with 3% BSA in PBS/Tween (100µl /well) for 1 hour at room temperature.
4. Wash the plates 3 times PBS/Tween as before and pat dry as before.
5. Apply (primary antibody) protein-specific antiserum (50µl/well) diluted from 1/50 in a doubling dilution series in PBS/Tween and incubate at room temperature for 90 minutes.
- 20 6. Wash plates as before (3 times quickly), followed up by 2 X 3 minute soaks (in PBS/Tween)
7. Apply diluted secondary antibody alkaline phosphatase conjugate. For anti-mouse Total IgG alkaline phosphatase conjugate (Goat Anti-Mouse IgG-AP, Southern Biotechnology Associates, Birmingham, AL. Cat. No. 1030-04) dilute 1/3000 in PBS/Tween and apply 50 µl per well and incubate at room temperature for 90 minutes.
- 25 8. Wash plates as in step 6.

9. Apply substrate. Dissolve one 5mg tablet of nitrophenyl phosphate (Sigma: kept in freezer) in 5ml of diethanolamine buffer. Apply 100  $\mu$ l per well. Cover with foil (a light-sensitive reaction) and leave at room temperature for 30 minutes. Read Optical densities (OD) at a wavelength of 405nm.
- 5 10. Plot curves of OD Vs dilution (log scale). Calculate end-point titres as the dilution giving the same OD as the mean of the OD obtained from the wells containing the 1/50 dilution of pre-immune serum.

10

15

### ELISA Plate format

2°	1/50	1/100	1/200	1/400	1/800	1/1600	1/3200	1/6400	1/12800	1/25600	1/51200
1°	Duplicate										
Pre											
Pre											
Pre											

### Table Summary

- 5     **Pre**     Replicate wells of pooled pre-inoculation serum (50µl per well) diluted to 1/50 are included on every plate in order for end point titres to be calculated.
- 2°     Is a blank control well to which no secondary antibody conjugate is applied. PBS/Tween by itself is applied instead
- 1°     Is a blank control well to which no primary antibody is applied. PBS/Tween by itself is applied instead
- 10     **Duplicate**     Each serum is analysed in duplicate
- The dilution series used is indicated (see first row of table). Beginning with a 1/50 dilution, sera are diluted two-fold in PBS/Tween in doubling dilution series as indicated.

15

### Protein Immunisation data

#### ID-65 and ID-83

The ID-65 and ID-83 vaccines were composed of the target proteins in phosphate buffered saline/10% glycerol mixed with aluminium hydroxide (alum)

(Imject® Alum, Pierce, Rockford, Ill.). Each dose of vaccine contained 20 µg of purified protein in 100 µl of PBS/10% glycerol, mixed with 50 µl of alum. A group of 6-8 week old CBA/ca mice (Harlan, UK) were immunised subcutaneously with the ID-65 and ID-83 vaccine and again 4 weeks later. A control group received a  
 5 150 µl dose of PBS/10% glycerol (2:1) with alum. All groups consisted of 6 mice. Mice were tail bled at 5 weeks post primary vaccination to obtain sera. The presence of total Immunoglobulin G (IgG) antibodies to the ID-65 and ID-83 protein in sera was determined by enzyme-linked immunosorbent assay (ELISA), using the purified protein as the coating antigen. ELISA was also performed using sera obtained at 6  
 10 weeks post-primary vaccination from the PBS/10% glycerol immunised control group.

NOTE: ELISA plates were coated with the ID-65 or ID-83 proteins at a concentration of 1 µg/ml.

15

#### **Protein Vaccination -ELISA results for ID-65 and ID-83**

Mice (6 per group) were immunised with two doses of the ID-65 and ID-83 vaccines with a four week interval. Mice were tail bled at 5 weeks post primary  
 20 vaccination to obtain sera. The Immunoglobulin G (IgG) titres to the vaccine protein component in sera were determined by enzyme-linked immunosorbent assay (ELISA), using the purified ID-65 and ID-83 proteins as the coating antigen. Subsequent to optimisation, ELISA plates were coated at a concentration 1µg/ml for both the purified ID-65 and ID-93 proteins. Total IgG titres were measured against  
 25 pre-immune serum (1/50 dilution). The results are shown in Table 2 and graphically in Figure 8.

30

Table 2

Serum (Group)	ID-65 + Alum (n=6)	PBS + Alum (n=6)	ID-83 + Alum (n=6)	PBS + Alum (n=6)
Coating antigen	ID-65		ID-83	
Bleed	5 weeks	5 weeks	5 weeks	5 weeks
Total IgG Titres (mouse 1-6)	7535763	965	82081	61
	1557649	90	50027	50
	3319737	108	154670	80
	1832259	176	57901	96
	8794360	371	66497	125
	1445728	0	49928	0
Average	4080916	285	76851	69
Standard Deviation	3258818	355	39985	43

5

# Protein Immunisation and Challenge data (ID-93)

## ID-93

10 The ID-93 vaccine was composed of the target ID-93 protein in phosphate buffered saline/10% glycerol mixed with aluminium hydroxide (alum) (Imject® Alum, Pierce, Rockford, Ill.). Each dose of vaccine contained 25 µg of purified protein in 100 µl of PBS/10% glycerol, mixed with 100 µl of alum. A group of 6-8 week old CBA/ca mice (Harlan, UK) were immunised subcutaneously with the ID-93 vaccine and again 4 weeks later. A control group received PBS/10% glycerol with alum. Both groups consisted of 6 mice. Mice were challenged at 7 weeks (unless otherwise stated). Mice were injected intraperitoneally (i.p.) with 5 X 10<sup>6</sup> bacteria diluted in

15

0.5 ml Todd-Hewitt broth. The challenged mice were observed daily for signs of illness. Deaths were recorded daily for 7 days. Survival data are shown in Table 3 and graphically in Figure 9.

5 Mice were tail bled at 3 weeks and 6 weeks post primary vaccination to obtain sera. The presence of total Immunoglobulin G (IgG) antibodies to the ID-93 protein in sera was determined by enzyme-linked immunosorbent assay (ELISA), using the pure ID-93 protein as the coating antigen. ELISA was also performed using sera  
10 obtained at 6 weeks post-primary vaccination from the PBS/10% glycerol immunised control group.

Note: ELISA plates were coated with the ID-93 protein at a concentration of 1  $\mu\text{g/ml}$ .

15 Table 3  
**ID-93 protein immunisation and GBS challenge experiment**

**Statistical analysis of Survival Times**

Group	PBS+ Alum	ID-93+ Alum
<u>Survival Times (hours)</u>	22.37	29.37
	22.37	35.12
	15.37	32.62
	28.03	32.62
	29.53	37.12
	26.53	27.87
<b>Mean</b>	24.03	32.45
<b>sd</b>	5.16	3.45
<b>p value</b>		0.01

20 **p value** refers to statistical significance when compared to unvaccinated controls.

# **Comment**

## **5 ID-93 (RS-70)**

Mice immunised with the ID-93-Alum vaccine exhibited significantly longer survival times when compared with the PBS-Alum control group.

10 (Statistical Significance was determined by the Mann-Whitney U test using a 95 % confidence level ( $p > 0.05$ ).

## **Protein Vaccination -ELISA results for ID-93**

15 Mice (6 per group) were immunised with two doses of the ID-93 vaccine with a four week interval. Mice were tail bled at 3 weeks and 6 weeks post primary vaccination to obtain sera. The Immunoglobulin G (IgG) titres to the vaccine protein component in sera were determined by enzyme-linked immunosorbent assay (ELISA), using the purified ID-93 protein as the coating antigen. Subsequent to optimisation, ELISA plates were coated with the purified ID-93 protein at a concentration of 1  $\mu\text{g/ml}$ .  
20 Total IgG titres were measured against pre-immune serum (1/50 dilution). The results are shown in Table 4 and graphically in Figure 10.

Table 4

Serum Group	<u>ID-93 + Alum(n=6)</u>		<u>PBS/10% glycerol (n=6)</u> (control)	
Coating antigen	<u>ID-93</u>	<u>ID-93</u>	<u>ID-93</u>	<u>ID-93</u>
Bleed	<u>3 weeks</u>	<u>6 weeks</u>	<u>3 weeks</u>	<u>6 weeks</u>
<u>Total IgG Titres</u> (mouse 1-6)	87196	3000000	39	100
	99544	8000000	31	16
	19620	2000000	31	79
	34724	10000000	59	48
	59990	10000000	24	328
	30041	4000000	13	40
<u>Average</u>	55186	6166667	33	102
<u>Standard error</u>	32654	3600926	15	115

5

#### Protein Immunisation data ID-89 and ID-96

10 The ID-89 and ID-96 vaccines were composed of the target proteins in phosphate buffered saline/10% glycerol mixed with TitreMax Gold adjuvant (Sigma, Missouri, USA) according to the manufacturers instructions. The ID-89 vaccine contained 25  $\mu\text{g}$  of purified protein 50  $\mu\text{l}$  of PBS/10% glycerol, mixed with 50  $\mu\text{l}$  of TitreMax Gold. The ID-96 vaccine contained 12.5  $\mu\text{g}$  of purified protein 50  $\mu\text{l}$  of PBS/10% glycerol, mixed with 50  $\mu\text{l}$  of TitreMax Gold. Groups of 6-8 week old CBA/ca mice

15 (Harlan, UK) were immunised subcutaneously with the ID-89 and ID-96 vaccines and again 4 weeks later. A control group received a 100  $\mu\text{l}$  dose PBS/10% glycerol with TitreMax Gold (1:1). Both groups consisted of 6 mice. Mice were tail bled at 3



weeks and 6 weeks post primary vaccination to obtain sera. The presence of total Immunoglobulin G (IgG) antibodies to the ID-65 and ID-83 protein in sera was determined by enzyme-linked immunosorbent assay (ELISA), using the purified protein as the coating antigen. ELISA was also performed using sera obtained at 3 weeks and 6 weeks post-primary vaccination from the PBS/10% glycerol immunised control group.

Note: ELISA plates were coated with the ID-89 or ID-96 proteins at a concentration of 1  $\mu$ g/ml and 3  $\mu$ g/ml respectively.

#### **Protein Vaccination -ELISA results for ID-89 and ID-96**

Mice (6 per group) were immunised with two doses of the ID-89 and ID-96 vaccines with a four week interval. Mice were tail bled at 3 weeks and 6 weeks post primary vaccination to obtain sera. The Immunoglobulin G (IgG) titres to the vaccine protein component in sera were determined by enzyme-linked immunosorbent assay (ELISA), using the purified ID-65 and ID-83 proteins as the coating antigen. Subsequent to optimisation, ELISA plates were coated with purified ID-89 and ID-96 protein at a concentration 1  $\mu$ g/ml and 3  $\mu$ g/ml respectively. Total IgG titres were measured against pre-immune serum (1/50 dilution). ELISA was also performed on both proteins using sera obtained at 3 weeks and 6 weeks post-primary vaccination from the PBS/10% glycerol immunised control group. Results are shown in tables 5a and 5b and graphically in Figure 11.

Table 5a

Serum	ID-89 + TitreMax Gold (n=6)		ID-96 + TitreMax Gold(n=6)	
Coating antigen	ID-89		ID-96	
Bleed	<u>3 weeks</u>	<u>6 weeks</u>	<u>3 weeks</u>	<u>6 weeks</u>
<u>Total IgG</u> <u>Titres</u> <u>(mouse 1-</u> <u>6)</u>	146940	1000000	190371	10000000
	89672	1000000	212505	10000000
	173532	2000000	167613	5000000
	85161	751210	110378	5000000
	88956	551281	142614	1000000
	27880	2000000	191085	1000000
Average	102024	1217082	169094	5333333
Standard Deviation	51451	629364	37341	4033196

Table 5b

Serum	<u>PBS/10%glycerol (n=6)</u>		<u>PBS/10%glycerol (n=6)</u>	
Coating protein	<u>ID-89</u>		<u>ID-96</u>	
Bleed	<u>3 weeks</u>	<u>6 weeks</u>	<u>3 weeks</u>	<u>6 weeks</u>
<u>Total IgG Titres (mouse 1-6)</u>	3	7	33	31
	8	18	77	62
	29	31	77	1
	34	4	52	29
	0	2	125	31
	5	1	113	0
Average	13	11	80	26
Standard deviation	15	12	35	23

#### **Example 4**

##### **5 Conservation and variability of candidate vaccine antigen genes among different isolates of Group B Streptococci**

An initial Southern blot analysis was carried out to determine cross-serotype conservation of novel Group B Streptococcal genes isolated using the LEEP system unless stated otherwise. Analysing the serotype distribution of a target gene will also  
 10 determine their potential use as antigen components in a GBS vaccine. The Group B Streptococcal strains whose DNA was analysed as part of this study are listed in APPENDIX III

##### **15 Amplification and labelling of specific target genes as DNA probes for Southern blot analysis.**

Oligonucleotide primers were designed for each individual gene of interest derived using the LEEP system unless stated otherwise. The same primers already described in APPENDIX II were used to amplify corresponding gene-specific DNA probes. Specific gene targets were amplified by PCR using Vent DNA polymerase (NEB)  
 20 according to the manufacturers instructions. Typical reactions were carried out in a 100  $\mu$ l volume containing 50 ng of GBS template DNA, a one tenth volume of enzyme reaction buffer, 1  $\mu$ M of each primer, 250  $\mu$ M of each dNTP and 2 units of Vent DNA polymerase. A typical reaction contained an initial 2 minute denaturation at 95°C, followed by 35 cycles of denaturation at 95°C for 30 seconds, annealing at  
 25 the appropriate melting temperature for 30 seconds, and extension at 72°C for 1 minute (1 minute per kilobase of DNA being amplified). The annealing temperature was determined by the lower melting temperature of the two oligonucleotide primers. The reaction was concluded with a final extension period of 10 minutes at 72°C.

All PCR amplified products were extracted once with phenol chloroform (2:1:1) and once with chloroform (1:1) and ethanol precipitated. Specific DNA fragments were isolated from agarose gels using the QIAquick Gel Extraction Kit (Qiagen). For use as DNA probes, purified amplified gene DNA fragments were labelled with digoxigenin using the DIG Nucleic Acid Labelling Kit (Boehringer Mannheim) according to the manufacturer's instructions.

#### **Southern blot hybridisation analysis of Group B Streptococcal genomic DNA**

Genomic DNA had previously been isolated from all strains of Group B Streptococci which were investigated for conservation of LEEP-derived (unless stated otherwise) gene targets. Appropriate DNA concentrations were digested using either *Hin* DIII or *Eco* RI restriction enzymes (NEB) according to manufacturer instructions and analysed by agarose gel electrophoresis. Following agarose gel electrophoresis of DNA samples, the gel was denatured in 0.25M HCl for 20 minutes and DNA was transferred onto Hybond™ N<sup>+</sup> membrane (Amersham) by overnight capillary blotting. The method is essentially as described in Sambrook *et al.* (1989) using Whatman 3MM wicks on a platform over a reservoir of 0.4M NaOH. After transfer, the filter was washed briefly in 2x SSC and stored at 4°C in Saran wrap (Dow chemical company).

Filters were prehybridised, hybridised with the digoxigenin labelled DNA probes and washed using conditions recommended by Boehringer Mannheim when using their DIG Nucleic Acid Detection Kit. Filters were prehybridised at 68°C for one hour in hybridisation buffer (1% w/v supplied blocking reagent, 5x SSC, 0.1% v/v N-lauryl sarcosine, 0.02% v/v sodium dodecyl sulphate[SDS]). The digoxigenin labelled DNA probe was denatured at 99.9°C for 10 minutes before being added to the hybridisation buffer. Hybridisation was allowed to proceed overnight in a rotating Hybaid tube in a Hybaid Mini-hybridisation oven. Unbound probe was removed by washing the filter twice with 2x SSC- 0.1% SDS for 5 minutes at room

temperature. For increased stringency filters were then washed twice with 0.1x SSC-0.1% SDS for 15 minutes at 68°C. The DIG Nucleic Acid Detection Kit (Boehringer Mannheim) was used to immunologically detect specifically bound digoxigenin labelled DNA probes.

### Results of Southern blot analysis

Unless otherwise stated, all genomic digests and their corresponding Southern blots followed an identical lane order as described in Table 6 below.

**Table 6**

Lane	1	2	3	4	5	6	7
Strain	1 kb molecular	515	A909	SB35	H36B	18RS21	1954/92
Serotype	Weight Marker	Ia	Ia	Ib	Ib	II	II

Lane	8	9	10	11	12	13	14
Strain	118/158	97/0057	BS30	M781	97/0099	3139	1169-NT
Serotype	II	II	III	III	III	IV	V

Lane	15	16	17	18	19	20
Strain	GBS 6	7271	JM9	Group A Streptococcus	<i>Streptococcus pneumoniae</i>	1 kb molecular
Serotype	VI	VII	VIII	-	14	Weight Marker

For comparative purposes, it was decided to analyse the serotype distribution of the GBS *rib* gene, which encodes the known protective immunogen Rib. Rib has previously been shown to be present in serotype III and some strains of serotype II but not in serotypes Ia or Ib (Stalhammar-Carlén *et al.*, *J. Exp. Med.* 177: 1593-1603 (1993)).

Confirmation of this pattern would not only give increased confidence in interpreting subsequent results, it would also determine if a *rib* gene homologue was present in the remaining GBS serotypes being investigated here. Primers designed for the amplification of *rib* for use as a gene probe in Southern blot analysis are described in APPENDIX II.

Table 7 – Lane order for Figure 12 (*rib* gene Southern blot analysis)

Lane	1	2	3	4	5	6	7
Strain	1 kb molecular	515	A909	SB35	H36B	18RS21	1954/92
Serotype	Weight Marker	Ia	Ia	Ib	Ib	II	II

Lane	8	9	10	11	12	13	14
Strain	118/158	97/0057	BM110	BS30	M781	97/0099	3139
Serotype	II	II	III	III	III	III	IV

Lane	15	16	17	18	19	20
Strain	1169-NT	GBS 6	7271	JM9	Group A Streptococcus	<i>Streptococcus pneumoniae</i>
Serotype	V	VI	VII	VIII	—	14

#### Rib (Figure 12) Comment

5 The Southern blot analysis shown in Figure 12 indicates that the *rib* gene is not conserved across all GBS serotypes. *rib* appears to be absent from all serotype Ia and Ib strains (lanes 2 to 5) and from strains 118/158 and 97/0057 of serotype II (lanes 8 and 9). However, *rib* would appear to present in strains 18RS21 and 1954/92 of serotype II (lanes 6 and 7) and in all strains of serotype III (lanes 10 to 13). This is in agreement with previously published data (Stalhammar-Carlemalm *et al.*, 1993 [*supra*]). *rib* would also appear to be present in strains representing serotypes VII and VII (lanes 17 and 18) but was absent from strains representing serotypes IV, V and V (lanes 14 to 16) as well as the control strains (lanes 19 and 20). The *rib* gene probe did hybridise with lower intensity to genomic DNA fragments from strains representing serotypes Ia, Ib, IV, VI, VII and serotype II strains 118/158 and 97/0057. This may indicate the presence of a gene in these strains with a lower level of homology to *rib*. These hybridising DNA fragments may contain a homologue of the GBS *bca* gene encoding the Ca protein antigen which has been shown to be closely homologous to the Rib protein (Wastfelt *et al.*, 20 *J. Biol. Chem.* 271:18892-18897 (1996)). If this is the case, it would be in agreement with previous work which showed all strains of serotypes Ia, Ib, II and III to be positive for one the two proteins (Stalhammar-Carlemalm *et al.*, 1993 [*supra*]). However, the apparent variable distribution of the *rib* gene amongst different GBS



serotypes, makes it a less than ideal candidate for use in a GBS vaccine that is cross-protective against all serotypes.

5

#### **ID-65 (Figure 13) Comment**

The Southern blot analysis described in Figure 13 indicates that gene ID-65 is conserved across all GBS serotypes. The gene probe hybridised specifically to a *Hin* DIII-digested genomic DNA fragment of approximately 3.0 kb in DNA digests from all GBS representatives, and was absent from both the control strains (lanes 18 and 19). This would suggest that the ID-65 gene is conserved across all GBS serotypes (and strains) at both the gene and locus level. The ID-65 DNA probe also hybridised weakly to the 1.636 bp molecular weight marker (the 1 kb DNA ladder from NEB was used to estimate DNA fragment sizes in all Southern blot analyses).

15

#### **ID-89 (Figure 14) Comment**

The Southern blot analysis described in Figure 14 indicates that gene ID-89 may not be conserved across all GBS serotypes. A 4.0 kb *Hin*DIII-digested genomic DNA fragment from 12 out of 16 GBS strains hybridised specifically to the ID-89 gene probe. In addition, a 3.25 kb *Hin*DIII-digested genomic DNA fragment from the GBS strain Ib (SB35) [lane 4] also hybridised specifically with the ID-89 gene probe. However, the ID-89 gene probe did not hybridise to digested genomic DNA fragments from strains Ia (515) [lane 2], IV (3139) [lane 13] and V (1169-NT) [lane 14], suggesting that these strains do not possess a ID-89 gene homologue.

25

#### **ID-93 (Figure 15) Comment**

The Southern blot analysis described in Figure 15 indicates that gene ID-93 is conserved across all GBS serotypes. The gene probe hybridised specifically to a *Hin*

DIII-digested genomic DNA fragment of approximately 3.25 kb in DNA digests from all GBS representatives, and was absent from both the control strains (lanes 18 and 19). This would suggest that the ID-93 gene is conserved across all GBS serotypes (and strains) at both the gene and locus level.

5

**ID-96 (Figure 16) Comment**

The Southern blot analysis described in Figure 16 indicates that gene ID-96 is conserved across all GBS serotypes. The gene probe hybridised specifically to a *Eco* RI-digested genomic DNA fragment of approximately 12.0 kb in DNA digests from all GBS representatives, and was absent from both the control strains (lanes 18 and 19). This would suggest that the ID-96 gene is conserved across all GBS serotypes (and strains) at both the gene and locus level.

10

## APPENDIX I

### ID-65

Forward Primer

5' - cggatccgccaccatgGCGGATCAA $\Delta$ CTACATCGGTTC - 3'

Reverse Primer

5' - ttgcggccgcGTTGGGATAACTAGTCGGTTTAGTCG

Length (including restriction sites) = 1541bp

10 Incorporating 1515bp of gene-specific sequence encoding 505 amino acids of the putative mature protein.

Annealing temperature for PCR amplification = 60°C

Sequence predicted to encode a signal peptide was omitted from amplified product

### 15 ID-66

Forward Primer

5' - cggatccgccaccatgAATCTTTATTTCCATAGTACTCCCTTGC - 3'

Reverse Primer

5' - ttgcggccgcAAAATGATCAGTTTGAGGGTAAAAGAG - 3'

20

Length (including restriction sites) = 767bp

Incorporating 747bp of gene-specific sequence encoding 247 amino acids of the putative mature protein.

Annealing temperature for PCR amplification = 60°C

25 Sequence predicted to encode a signal peptide was omitted from amplified product

## APPENDIX II

**ID-65**

Forward Primer

5' - catgccatgGCGGATCAAACACTACATCGGTTC - 3'

5 Reverse Primer

5' - ttgcggccgcGTTGGGATAACTAGTCGGTTTAGTCG

Length (including restriction sites) = 1534bp

10 Incorporating 1515bp of gene-specific sequence encoding 505 amino acids of the putative mature protein.

Annealing temperature for PCR amplification = 60°C

**ID-83**

15 Forward Primer

5' - catgccatggcaAAAATAGTAGTACCAGTAATGCCTC - 3'

ReversePrimer

5' - ttgcggccgcCTCTGAAATAGTAATTTGTCCG - 3'

20 Length (including restriction sites) = 626bp

Incorporating 624bp of gene-specific sequence encoding 208 amino acids of the putative mature protein.

Annealing temperature for PCR amplification = 52°C

25

**ID-89**

Forward Primer

5' - catgccatgggaAAGAAAGCAAATAATGTCAGTCC - 3'

**Reverse Primer**

5' – ttcgggccgcATTGGGTGTAAGCATTTTTTC –3'

Length (including restriction sites) = 990bp

5 Incorporating 969bp of gene-specific sequence encoding 323 amino acids of the putative mature protein.

Annealing temperature for PCR amplification = 54°C

**ID-93**

**Forward Primer**

10 5' – catgccatgggaACTGAGAACTGGTTACATACTAAAG – 3'

**ReversePrimer**

5' – ttcgggccgcATTAGCTTTTTCAACAATTTCTC – 3'

Length (including restriction sites) = 759bp

15 Incorporating 744bp of gene-specific sequence encoding 248 amino acids of the putative mature protein.

Annealing temperature for PCR amplification = 51°C

**ID-96**

20 **Forward Primer**

5' – ctagtagccgATGTTTTCGTGGGAAAG – 3'

**ReversePrimer**

5' – ttcgggccgcATAAGATTTAACAATACCAAGTAATATAGC – 3'

Length (including restriction sites) = 944bp

25 Incorporating 921bp of gene-specific sequence encoding 307 amino acids of the putative mature protein.

Annealing temperature for PCR amplification = 53°C

*rib* (control)

Forward primer

5' - ggggtaccggccaccATGGCTGAAGTAATTTTCAGGAAGT -3'

5 Reverse primer

5' - cggaattccgTTAATCCTCTTTTTTCTTAGAAACAGAT

Length (including restriction sites) = 3559bp

10 Incorporating 3531bp of gene-specific sequence encoding 1177 amino acids of the mature protein.

Annealing temperature for PCR amplification = 55°C

**APPENDIX III**

15 Listed below are the details (serotype and strain designation) of Group B Streptococcus strains whose DNA was analysed for gene conservation

		SEROTYPE	STRAIN
20	Ia	515	
	Ia	A909	
	Ib	SB35	
	Ib	H36B	
	II	18RS21	
25	II	1954/92	
	II	118/158	
	II	97/0057	
	III	BM110	
	III	BS30	
30	III	M781	
	III	97/0099	
	IV	3139	

A group A Streptococcal strain (serotype M1, strain NCTC8198) and *Streptococcus pneumoniae* (serotype 14) were also included in the analysis for control purposes.